

BLACK GLASS STANDARDS



Low reflectance glass standard

Categories: PHYSICAL STANDARDS, REFLECTANCE STANDARDS



- Product Information**

Highly polished black glass is useful as a standard for gloss, low reflectance in directional/directional geometry, and also establishing the efficiency of the 'sphere zero' in specular excluded geometry using an integrating sphere.

Options/Specifications

Calibration Calibrated 360-830 nm 8°/hemispherical Geometry, or uncalibrated

Dimensions 2"/52 mm square; custom sizes available

- Ordering**

BG-02 *2" square black glass standard with holder*

Custom *Custom products are our specialty!
Contact us for more information.*



Avian Technologies

Solutions for Optical Spectroscopy Materials, Coatings, and Standards

News

NEWS

LC Instru Represents Avian Technologies

Avian Technologies is pleased to welcome our newest representative, **LC Instru**. If you are interested in purchasing Avian Technologies products in France, give them a call!

Visit our **international sales** page for more information on purchasing Avian Technologies' products if you are located outside the US.

Contact :



LC-Instru SAS / JASCO France SAS

ZI de l'Eglantier

17 rue des Cerisiers

91090 LISSES

Tél : 09 72 50 12 35 ou 01 64 97 09 60

info@lc-instru.fr ou jascofrance@jascofrance.fr