


REFLECTANCE MINI-ARRAY & MICRO-ARRAY STANDARDS



Small, portable arrays of gray scale standards



Categories: [PHYSICAL STANDARDS](#), [REFLECTANCE STANDARDS](#)



- [Product Information](#)

Avian Technologies offers micro-arrays of our popular [grey scale standards kits](#). These arrays are ideal for use with portable spectrophotometers for the Vis-NIR range, which are typically able to measure a very small area. Our micro-arrays are more convenient than large sets of standards, making them more suitable for field work.

Arrays include either 4 or 8 standards. The standards are bulk calibrated, meaning that larger standards are calibrated and then cut into small parts and assembled as an array. Measurement data is provided at 1 and 20 nm increments.



Product Options/Specifications

Individual Standard Dimensions	12.5 mm x 12.5 mm
Array Dimensions	30 mm x 60 mm
Standards Available	2%, 5%, 10%, 20%, 40%, 60%, 80%, 99% nominal reflectance at 600 nm

- Ordering**

MNA-FSS08-C *Set of 8 standards (2, 5, 10, 20, 40, 60, 80, 99%) in one holder*

MNP-FSS04-C *Set of 4 standards (5, 10, 40, 99%) mounted in 60 mm x 30 mm holder*

MNP-FSS08-C *Set of 8 standards (2, 5, 10, 20, 40, 60, 80, 99%) mounted in 60 mm x 30 mm holder*

Custom *Custom products are our specialty!
Contact us for more information.*



Avian Technologies

Solutions for Optical Spectroscopy Materials, Coatings, and Standards

News

NEWS

LC Instru Represents Avian Technologies

Avian Technologies is pleased to welcome our newest representative, **LC Instru**. If you are interested in purchasing Avian Technologies products in France, give them a call!

Visit our **international sales** page for more information on purchasing Avian Technologies' products if you are located outside the US.

Contact :



LC-Instru SAS / JASCO France SAS

ZI de l'Eglantier

17 rue des Cerisiers

91090 LISSES

Tél : 09 72 50 12 35 ou 01 64 97 09 60

info@lc-instru.fr ou jascofrance@jascofrance.fr