



---

Nominal reflectances at 600 nm 2%, 5%, 10%, 20%, 40%, 60%, 80%, 99%\*

---

Calibration Calibrated or uncalibrated

\*Tolerance is +/- 1% at 600 nm for 2%, 5%, and 99% standards. All others have a tolerance of +/- 2%.

- **Ordering**

<b>EFGS-xx-01</b>	<i>1.25" diameter Fluorilon Grey Scale, Encapsulated Available calibrated or uncalibrated "xx" = nominal reflectance</i>
<b>EFGS-xx-02</b>	<i>2.00" diameter Fluorilon Grey Scale, Encapsulated Available calibrated or uncalibrated "xx" = nominal reflectance</i>
<b>EFSS-04-01</b>	<i>1.25" diameter, Set of 4 Encapsulated Grey Scale Reflectance Standards: 2%, 99%, plus two of your choice Available calibrated or uncalibrated, traceable to NIST/NRC, boxed</i>
<b>EFSS-04-02</b>	<i>2.00" diameter, Set of 4 Encapsulated Grey Scale Reflectance Standards: 2%, 99%, plus two of your choice Available calibrated or uncalibrated, traceable to NIST/NRC, boxed</i>
<b>EFSS-08-01</b>	<i>1.25" diameter, Set of 8 Encapsulated Grey Scale Reflectance Standards: 2%, 5%, 10%, 20%, 40%, 60%, 80%, 99% Available calibrated or uncalibrated, traceable to NIST/NRC, boxed</i>
<b>EFSS-08-02</b>	<i>2.00" diameter, Set of 8 Encapsulated Grey Scale Reflectance Standards: 2%, 5%, 10%, 20%, 40%, 60%, 80%, 99% Available calibrated or uncalibrated, traceable to NIST/NRC, boxed</i>
<b>Custom</b>	<i>Custom products are our specialty! Contact us for more information.</i>



## Avian Technologies

Solutions for Optical Spectroscopy Materials, Coatings, and Standards

# News

### NEWS

## LC Instru Represents Avian Technologies

Avian Technologies is pleased to welcome our newest representative, **LC Instru**. If you are interested in purchasing Avian Technologies products in France, give them a call!

Visit our **international sales** page for more information on purchasing Avian Technologies' products if you are located outside the US.

### Contact :



LC-Instru SAS / JASCO France SAS

ZI de l'Eglantier

17 rue des Cerisiers

91090 LISSES

Tél : 09 72 50 12 35 ou 01 64 97 09 60

[info@lc-instru.fr](mailto:info@lc-instru.fr) ou [jascofrance@jascofrance.fr](mailto:jascofrance@jascofrance.fr)