


SPECULAR STANDARDS

First surface aluminum mirror standards

Categories: PHYSICAL STANDARDS, REFLECTANCE STANDARDS




- **Product Information**



Avian Technologies offers first surface specular standards with calibrations traceable to NIST/NRC. The mirrors are first surface aluminum with a magnesium fluoride overcoat that have been aged to enhance their long term stability.

Standards come in a protective holder.

Options/Specifications



Hemispherical reflectance factor - Avian Technologies Front Surface Specular Standard (click for full size)

If you require measurements outside of the standard calibration data, we will be happy to refer you to one of our associates to provide additional measurements.

Ordering

- SPFS-01d** *1.25" diameter front surface aluminum mirror, calibrated, in holder*
- SPFS-02s** *2" square front surface aluminum mirror, calibrated*
- Custom** *Custom products are our specialty!
Contact us for more information.*



Avian Technologies

Solutions for Optical Spectroscopy Materials, Coatings, and Standards

News

NEWS

LC Instru Represents Avian Technologies

Avian Technologies is pleased to welcome our newest representative, **LC Instru**. If you are interested in purchasing Avian Technologies products in France, give them a call!

Visit our [international sales](#) page for more information on purchasing Avian Technologies' products if you are located outside the US.

Contact :



LC-Instru

LC-Instru SAS / JASCO France SAS

ZI de l'Eglantier

17 rue des Cerisiers

91090 LISSES

Tél : 09 72 50 12 35 ou 01 64 97 09 60

info@lc-instru.fr ou jascofrance@jascofrance.fr